



FloorBed[®] 1

By Accora

Specification Sheet



Reduce the risk of injury in the event of a patient fall

Reduce the risk of caregiver injury

Enable independence and mobility

OVERVIEW

Order code:

NSB-0-FL1-200

SPECIFICATION

Overall dimensions:	925 mm x 2255 mm (36.4 x 88.8 in)
Mattress size (standard configuration):	900 mm x 2000 mm (35.4 x 78.7 in)
Mattress size (extended configuration):	900 mm x 2200 mm (35.4 x 86.6 in)
Safe working load:	185 kg (29 st)
Patient weight:	150 kg (24 st)
Lift range (excl. mattress):	71 mm - 650 mm (2.8 - 25.6 in)
Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg

TECHNICAL INFORMATION

Head and foot bumpers BMHNE-0-FL1-100	<ul style="list-style-type: none"> Consists of both foot and head bumper. Attachable to the frame by means of nylon straps and plastic buckle clips. Attachable to FloorBed® Safety Sleeve by means of a zip.
Standard bed lever STLEV-0-FL1-100	<ul style="list-style-type: none"> Can be placed in several different positions on the bed frame. Provides support for bed mobility and transfers.
Rotating bed lever RTLEV-0-FL1-000	<ul style="list-style-type: none"> The additional length provides extra support with lying to sitting movement. When rotated, this provides additional support from sit to stand, and may enable to "pull" themselves up rather than push.
Integrated padded side rails – standard length SDR-0-FL1-000	<ul style="list-style-type: none"> Full length support. Easy lock and lower mechanism that locks securely in position. Side rails do not impede on transfers as they are below the mattress when lowered. 2 in 1 bumper and side rail. Compatible with headboard / footboard bumpers.
Integrated padded side rails – for use when FloorBed® is extended by 20 cm SDREX-0-FL1-000	<ul style="list-style-type: none"> Total length 220 cm.
FloorBed® safety sleeve SLE-0-FL1-000	<ul style="list-style-type: none"> Safety sleeve may reduce risk of entrapment by reducing gaps when raising the knee break and back rest functions (protecting the platform mechanism). Risk of entrapment is further reduced when connected to the footboard and headboard bumpers.
Support bar / Self assist pole LIFOL-0-FL1-000	<ul style="list-style-type: none"> Support for sitting up or changing position in bed.
High safety mat SAFTH-0-FL1-000	<ul style="list-style-type: none"> 22 cm high safety mat - top of safety mat is same height as the top of the mattress when a 15 cm mattress is used, and FloorBed® is at its lowest position. Folds for storage.
High safety mat with integrating sliding sheets SAFTHF-0-FL1-000	<ul style="list-style-type: none"> 22 cm high safety mat - top of safety mat is same height as the top of the mattress when a 15 cm mattress is used, and FloorBed® is at its lowest position. Integrated slide sheets for ease of transfer. Folds for storage.
Multicare 3000 Visco Mattress VISM-0-FM1-000	<ul style="list-style-type: none"> Layered soft memory foam technology for improved pressure reduction. Base layer of high density foam for support. Robust foam structure for longevity.
Mattress infill MTINF-0-FL1-000	<ul style="list-style-type: none"> Infills to fit between mattress and side rails.